Type of material	Destination	Size	Language	Variant	Date
PIL	ABIOJOINT EVO		--	ABIOGEN	21/04/26
Dimensions (mm)	N° user	Color 1	Color 2	Color 3	Color 4
500x600	2	P 1797	BLACK	P COOL GRAY	
Packaging/Formato	IMPIANTI DI PROPRIETÀ	ABIOGEN Pharma S.p.A.			VIETATA LA MANOMISSIONE RENDERE DOPO LA STAMPA
					IT11021-05

INSTRUCTIONS FOR USE

NAME
Abiojoint EVO

DESCRIPTION
Abiojoint EVO is a sterile, biodegradable, transparent gel of non-animal cross linked hyaluronic acid (HA synthesized by bio-fermentation) based on HA crosslinking technology EVENOR TOUCH®. The device is supplied with a prefilled glass syringe filled with 2 ml of gel. The content of the syringe is sterilized using moist heat. The syringe is packed individually in a blister. The product is for single use only.

PERFORMANCE AND CHARACTERISTICS
Abiojoint EVO is a cross-linked HA gel for intra-articular use that achieves its effect through viscosupplementation. Viscosupplementation is a process whereby the physiological and rheological states of the impaired joint are restored. Viscosupplementation is achieved by injection of exogenous hyaluronic acid into diarthro-ideal joints, with the aim of reducing the articular pain in the treated knee and improving the knee mobility. Abiojoint EVO acts temporarily only at the joint where it is injected, without exerting any systemic action. The intended clinical benefits of the cross-linked HA gel for intra-articular use are:

- Restoration of the correct lubrication of the joint;
- Reduction of pain, stiffness and loss of flexibility thanks to the shock absorber and cushioning properties typical of the physiological synovial fluid.

Abiojoint EVO is characterized by a neutral pH (6.5 – 7.5) and by an Osmolarity between 270 and 380 mOsm/L, which is suitable for an intra-articular injection. Hyaluronate is a component of synovial fluid, which is responsible for its elastoviscosity. The mechanical (elastoviscous) properties of Abiojoint EVO are similar to those of synovial fluid and superior to those of unmodified hyaluronate solutions of comparable concentration.

Abiojoint EVO has an elastosity (storage modulus G') at 1 Hz between 15 and 60 Pascals (Pa) and a phase angle (tan δ) between 0.2 and 0.8. Elastic and viscous moduli of knee synovial fluid of 18 to 27 years old humans measured with comparable method at 2.5 Hz are G' = 117 Pa and G'' = 45 Pa.

INDICATIONS
Abiojoint EVO is a temporary replacement and supplement of the synovial fluid which thanks to its viscoelastic properties, promotes the reduction of pain and improvement of knee mobility. Abiojoint EVO is intended only for intra-articular use by a physician. Abiojoint EVO is indicated for viscosupplementation in knee joint. Patient population are adults (> 18 years) male and female.

DOSAGE GUIDELINES
The recommended treatment regimen is intra-articular injection using either one injection (2ml) of Abiojoint EVO or two injections (2 ml each) of Abiojoint EVO one week apart. The results of clinical investigation showed that treatment with Abiojoint EVO was well tolerated in both regimens of administration (one or two I.A. injections). The use of the two-injection regimen was not associated with an increased risk of adverse events.

CONTRAINDICATIONS

- Do not use the same device after the first use.
- Do not use the same device with two different patients.
- Do not use in infected or severely inflamed joints or in patients having skin diseases or infections in the area of the injection site.
- Do not use if there is a large intra-articular effusion prior to the injection or if venous or lymphatic stasis is present in the relevant limb.
- Do not use if the patients have bleeding disorders.
- Do not use in patients with known allergy to sodium hyaluronate preparation.
- Do not use in pregnant or breastfeeding women.
- Do not use in children.

WARNINGS AND PRECAUTIONS

- Carefully read these instructions for use before using the product.
- Adhere to the indicated dosage guideline. An excess quantity of sodium hyaluronate should be avoided and the patient should be closely monitored, as the adverse events identified in the following section may occur.
- Do not concomitantly use disinfectants containing quaternary ammonium salts for skin preparation because hyaluronan can precipitate in their presence.
- The use of the operating procedure other than those published by the manufacturer, or the use of the accessory device not recommended by the manufacturer can result in patient injury or death.
- Physicians must discuss all potential risks of intra-articular injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.
- Due to the absence of clinical data in patients who had arthroplasty at the target knee, physicians must carefully evaluate the potential benefits and risks of intra-articular injection in this population.
- It is recommended to remove joint effusion, if present, before injecting the product.
- Abiojoint EVO must be administered by a physician only.
- Injection procedures are associated with a risk of infection. Aseptic technique and standard practice to prevent cross-infections are to be observed.
- Abiojoint EVO must be used in sterile operating conditions after adequate preparation of the site to be treated.
- It is important to use appropriate needle. Suitable sterile needles are specified in the following section "Instructions for Use".
- Do not use the product if the package is damaged or if the syringe or the blister lid are damaged or opened before use.
- Do not use the product if the contents of the syringe are cloudy.
- Abiojoint EVO must be used immediately after opening the blister lid.
- Do not perform intra-vascular injection or injection near to blood vessels.
- Abiojoint EVO is a single-use product and must not be re-sterilized.
- Abiojoint EVO must not be stored for later use once the blister lid has been opened as this might create a risk of infection for the patient. Reusing syringes already used to inject the product can result in the transmission of infectious agents, potentially endangering patients, physicians and staff.
- Particular attention must be paid in the connection, disconnection and disposal of the needle. Disposal of needles shall be done using the appropriate disposal containers suitable for sharps, according to national or local laws in order to avoid accidental puncture and any contamination risk.
- Disposal of empty syringes or containing any residual product shall be done using suitable medical waste containers according to national or local laws.

WARNINGS AND PRECAUTIONS FOR THE PATIENT

- Patients should carefully store the Implant Card associated with the device.
- As with any invasive joint procedure, it is recommended that the patient avoids any strenuous activities following the intra-articular injection, and resumes full activities within a few days.
- Transient pain or swelling may occur after the intra-articular injection.

ADVERSE EVENTS
As with any other intra-articular HA preparations, local pain, arthralgia, swelling, heat and redness may occur sporadically at the injection site. Such symptoms are generally mild and transient.

INSTRUCTIONS FOR USE

- Inject at room temperature.
- Strict adherence to aseptic techniques must be followed.
- The use of topical and subcutaneous anesthetic is not required but can be used at the discretion of the physician.
- Before Abiojoint EVO injection, remove synovial fluid or effusion, if present. Gently aspirate the joint.
- Do not use the same syringe for removing synovial fluid and for injecting Abiojoint EVO, but the same needle should be used.
- Remove the Abiojoint EVO syringe from its packaging by taking hold of the syringe barrel by the body without touching the plunger rod.
- Take particular care to remove the tip cap of the Abiojoint EVO syringe aseptically.
- Only 18 – 21 gauge needle is suitable for the injection of Abiojoint EVO. Secure the needle tightly while firmly holding the Luer hub to ensure a tight seal and prevent leakage during administration.
- Do not over tighten or apply excessive leverage when attaching the needle or removing the needle guard, as this may break the tip of the syringe.
- Inject Abiojoint EVO into the synovial space only.
- The physician should ensure proper penetration into joint synovial space prior to injecting Abiojoint EVO.
- Administer the whole volume of the syringe.
- Only the needles with standard luer lock connection assure a safe connection with syringe of Abiojoint EVO.
- When the treatment is completed, the healthcare institution or healthcare provider must fill out the Implant Card, as per instructions included on it, and provide it to the patient together with this instruction for Use.

CONTENTS

Composition of Abiojoint EVO:

Component	Quantity per syringe (2 ml)
Sodium hyaluronate crosslinked*	40 mg
Sodium chloride	18 mg
Potassium dihydrogen phosphate	0.04 mg
Disodium phosphate dihydrate	0.44 mg
Water for injection	Up to 2 mL

*Crosslinker = 1,4 Butanediol diglycidyl ether (BDDE)

NATURE OF THE CONTAINER

Abiojoint EVO is packaged in 2.25 ml Type I glass syringes filled with 2 mL of gel, closed with a tip cap on one side and with a plunger seal on the other side, both made of elastomeric materials. Each syringe is equipped with a luer lock, a plunger rod and a finger rest made of polypropylene.

DEVICE LIFETIME

The lifetime, meant like the performance of the product, is 3 months based on primary endpoint of clinical investigation. From the perspective of clinical effects duration, when Abiojoint EVO is injected, it replaces or integrates the natural HA; available data from clinical investigation demonstrate that positive effects, intended as significant reduction of the articular pain in the target knee, are up to 9 months.

Preclinical and literature data indicate that the persistence of the product, changes on the basis of disease severity, from few days (more severe inflammation) to 20-30 days. The device starts to be degraded in smaller fragments that still maintain viscosupplemental effects as they are in the range of 300 and 500 kDa. *In-vitro* tests show a degradation profile of Abiojoint EVO comparable to other HA-based preparation of similar marketed products.

SHELF LIFE, STORAGE AND HANDLING CONDITIONS

- Do not use after the expiry date indicated on the package.
- Store between 5°C and 25°C.
- Protect from freezing and sunlight.

DEFINITION OF SYMBOLS USED IN THE PRODUCT LABELLING

Symbol	Title of Symbol
	Manufacturer
	Date of manufacture
	Catalogue number
	Sterilized using steam or dry heat
	Sterile fluid path. Sterilized using steam or dry heat
	Batch code
	Medical Device
	Unique Device Identifier carrier
	Use-by date
	Do not resterilize
	Do not re-use
	Consult instructions for use
	Single sterile barrier system
	Sterilized using steam or dry heat
	Single sterile barrier system with protective packaging outside
	Do not use if package is damaged and consult instructions for use
	Fragile, handle with care
	Keep dry
	Temperature limit
	Keep away from sunlight

NOTICE TO THE USER/PATIENT

Any serious incident occurred in relation to the device should be reported to the manufacturer (reclamIMDEV@abioigen.it) and the relevant competent authority.

INFORMATION FOR THE USER/PATIENT

These instructions for use are available also to the following web address: <https://www.abioigen.it/it/medical-device/>

Summary of Safety and Clinical Performance (SSCP) is available in the European database (Eudamed), linked to the Basic UDI-DI. Website: <https://ec.europa.eu/tools/eudamed> Basic UDI-DI: 805651586P90040211001ZL

MANUFACTURER IDENTIFICATION

	ABIOGEN PHARMA SpA
	Via Meucci 36, 56121 Ospedaletto (PI) - ITALY
	www.abioigen.it

Instructions for Use version April 2026

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ABIOGEN PHARMA S.p.A.